



Handicare 2000

Twin rail curved stairlifts



The benefits of a stairlift from Handicare

If you have difficulty getting up and down stairs you are not alone. Every year thousands of people look for a solution to the problems posed by the difficulty they have using the staircase in their home.

Handicare's stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.


The decision to install a stairlift is sometimes a difficult one to take but when compared with the alternative of moving house, or living on one level of their house, many people have found it to be their preferred solution.

If you are thinking about installing a stairlift and your stairs turns corners, the Handicare 2000 offers you a safe and reliable solution. The stylish track system now gives you it's smoothest ever ride.

The Handicare 2000 stairlift suits most staircases and offers you all the basic functions that you would expect plus the option to upgrade to a powered folding footrest and/or powered seat swivel.

A tailor-made solution


Made to measure and perfectly finished it can be quickly installed at your home. After an informal visit we explain all the features of our stairlifts to you in detail.



Style Seat

Each Handicare stairlift is designed to suit your individual staircase by using advanced PhotoSurvey3D technology, to make quick and accurate measurements. We have all the options available to ensure the best possible fit for your staircase.

Your Handicare 2000 stairlift can be fitted with a choice of seats and our Handicare supplier will be able to advise you which solution is most suited to your needs.



Smart Seat



Simplicity Seat



Standard Features



Remote control:

Your Handicare 2000 stairlift comes with two remote controls as standard, so you can have one upstairs and one downstairs. This means you can send the lift up and down the track without having to sit on it. This can also help free up space in your hallway or landing.



Safety as standard:

Your Handicare 2000 stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstructions during travel. All Handicare stairlifts have a battery backup, which means it will continue to work in the event of a power cut.



Foldaway:

The stairlift folds away, leaving the staircase free for other users. The tailor made track is designed to maximise space on the stairway.





Easy to use:

The operating toggle is very easy to use, simply hold it towards the direction you wish to travel and your stairlift will move smoothly on its way.



Seat belt:

For added peace of mind, your Handicare 2000 is fitted with an easy to use retractable seat belt. If you prefer, you can choose the Velcro seat belt solution.



Diagnostic display:

Your Handicare 2000 stairlift has a clear diagnostic display, so that you can monitor the status of your stairlift, for example if a safety sensor has been activated.



Comfort features

Automatic hinge track

Some homes will have a doorway close to the bottom of the stairs. An automatically hinged track is available for these situations and means that the track will be folded away from the doorway to give clear access. When operating the stairlift from the toggle, the hinge will lift and lower automatically.



Powered seat swivel

When you arrive at the top of the stairs, you will need to swivel your seat round, so that you exit onto the landing safely. If you choose the powered option, the swivel is automatically operated by continuing to hold the toggle when you arrive at the top of the staircase.




Weight capacity

The standard maximum weight limit on the Handicare 2000 is 250 lbs. This weight limit can be increased to 300 lbs. by having the heavy duty kit fitted on your stairlift.

Powered footrest

The footrest and seat can be folded up to leave the staircase as clear as possible for other users to walk past. If you choose the powered footrest option it will remove the need for you to bend down to fold or unfold the footrest.



Stairlift colors

Fitted to the stairs and not the wall, it is unlikely you will need to redecorate following the installation of your Handicare 2000 stairlift. With a choice of upholstery and track colors, your stairlift will compliment most color schemes. For an additional fee, the rail can be finished in any color to suit your décor.

Perch seat

If you have a narrow staircase or if restricted movement in the knee or hip joint makes sitting painful, a perch seat may be the solution.

A perch seat means that the user is almost in a standing position when using the lift. There is a shortened seat and retractable seat belt which provide a stable support area for users to 'perch' on. The upholstery is available in a choice of six colors.




Handicare 2000 series: at a glance

Opzioni sedile	2000 Smart	2000 Simplicity	2000 Style	
Choice of Upholstery	✓	-	✓	A selection of six or three colors means you can choose upholstery to match your décor.
Powered Swivel	○	-	○	This makes turning the stairlift seat to get in and out on the landing as easy as possible.
Retractable Seat Belt	✓	✓	✓	This retractable belt can be clipped together on one side or centrally.
Velcro Seat Belt	○	○	○	This black belt has velcro on either end to hold it firmly in place.
Perch Seat	○	-	-	Ideal for users who have difficulty bending at the hip or knee joint.
Footrest Options	2000 Smart	2000 Simplicity	2000 Style	
Manual	✓	✓	✓	The footrest is folded and unfolded manually.
Powered from Seat	○	○	○	Folds the footrest away when the seat is folded up.
Powered from Arm	○	-	-	The footrest can be lifted and lowered by a switch under the arm whilst the user is seated.
Track Features	2000 Smart	2000 Simplicity	2000 Style	
Travels around corners and bends	✓	✓	✓	For a wide range of staircases, from ones that turn a corner to multi-flight or spiral.
Multiple charge points	○	○	○	Use your remote control to park on charge points.
Track available in any color	○	○	○	Customise your curved stairlift track to match your décor.
Zero intrusion into landing	○	○	○	For use when a doorway close to the top of the stairs would be obstructed by the track.
Powered hinge	○	○	○	Automatically moves the track while the stairlift moves leaving doorways unobstructed.
Weight capacity	2000 Smart	2000 Simplicity	2000 Style	
Standard Maximum Weight	250 lbs	250 lbs	250 lbs	
Max weight with Heavy Duty Kit	300 lbs	300 lbs	300 lbs	
Standards compliance	2000 Smart	2000 Simplicity	2000 Style	
EN 81-40: Stairlifts	✓	✓	✓	


✓ Available as standard ○ Optional - Not available

Smart Seat: Upholstery Color Choices



COCOA RUBY SAND SAPPHIRE SLATE JADE

Simplicity Seat: Upholstery Color




SAND

Style Seat: Upholstery Color Choices



CHERRY RED VANILLA BEIGE IVORY WHITE

Track Color Choices:



LIGHT GREY BROWN IVORY WHITE

Handicare also offers special orders so the track can be finished in any color to suit your décor.

Please be aware that actual colors may vary slightly from those shown.



Making everyday life easier

Handicare Vision: a 3D preview of a stairlift on your stairs

When it comes to stairlifts we want you to make the best choice possible. You can have a 3D preview of Handicare's stairlift range in the comfort of your own home, and see how each one would look in your hallway too, thanks to Handicare's Vision app. Vision is a special application we have developed to show you what a stairlift would look like on your stairs, and to help you make the right decision.



Quality and safety

All our new products are delivered with a comprehensive warranty.

A Handicare stairlift stands for high quality: quality in terms of safety, innovation, comfort and design. Furthermore, all our products are tested and approved by worldwide independent inspection services. This gives you the assurance that all our products meet the highest standards.



Experience

Convenience serves people... and lifting equipment of Handicare has been doing just that for more than 125 years.

With technical perfection, intelligent solutions and safe applications, Handicare's innovations have provided not only greater comfort, but also more attractive stairlifts... in your style, your taste.

Top quality in all respects make the Handicare stairlift a valuable asset for independence around the home. Many satisfied users experience that every day.



Handicare
2201 Hangar Place, Suite 200
Allentown, PA 18109
<http://handicarestairlifts-us.com>
customerservice.usa@handicare.com

Local 610-266-5260
Toll free 866-276-5438
Fax 610-266-5266