

Full control in your hands.

Android and Apple compatible.



Telephone

Stream phone calls directly from your Apple® smartphone (and with StreamLine Mic from your Android™ smartphone) to your hearing aids for complete convenience.



Music

Listen to your favorite music streamed directly into your hearing aids from your Apple device (and with StreamLine Mic from your Android device) for an outstanding sound experience.



TV

With StreamLine TV, audio from your television can be fed directly into your hearing aids in Dolby® Digital stereo sound. You can easily stream and adjust the volume through the myControl App.



Remote control

The myControl App lets you change hearing aid settings remotely and personalize your hearing aids to suit your individual preferences.

And if you don't have a smartphone, the miniPocket™ also provides highly convenient remote control.

● for Pure 312 ● for Pure 13 ● for Pure 10

Natural hearing at its best.



Pure 312 / Pure 13 / Pure 10



Speech in noise

Thanks to **Speech Quality**, Signia's hearing aids enable you to hear speech clearly in all situations through advanced noise suppression technologies. This makes it easier for you to understand your conversation partner, even in very noisy environments such as a restaurant or a conference.



Crystal clear sound

Sound Clarity offers our most natural, lifelike sound in all environments. This lets you enjoy crystal clear sound, even in the loudest environments.

Pure 312 / Pure 13



Hearing your own voice

Many hearing aid wearers do not like the sound of their own voice as it can sound loud and unnatural. As a world's first, Signia's hearing aids solve this: their **OVP™ (Own Voice Processing)** provides you with a natural sounding own voice combined with optimal hearing of all other sounds for our most pleasant, natural listening experience.



Conversation while moving

The **3D Classifier** provides advanced automatic control in all situations to let you enjoy our highest hearing performance. To enjoy optimal hearing even when moving, the hearing aid settings are automatically adjusted based on the **myControl App's motion data**.

Pure 312 Nx



Maximum discretion

Pure sophistication in a small package: with its small, elegant housing Signia's Pure® 312 Nx sits unobtrusively behind your ear. It can be remotely controlled via your smartphone for a maximum of convenience and discretion.



Fully featured

Packed with state-of-the-art technology, Pure 312 Nx offers our most natural sound and our highest hearing performance as well as superior connectivity.

IP68 rated** – resistant to moisture, sweat, and dirt.

Pure 13 Nx



Maximum energy efficiency

The reliable Pure® 13 Nx from Signia provides the longest audio streaming time in its class, so you can focus on the important things in life.



Fully featured

Packed with sophisticated, state-of-the-art technology, Signia's Pure 13 Nx offers our most natural sound and our highest hearing performance as well as superior connectivity.

IP68 rated** – resistant to moisture, sweat, and dirt.

Pure 10 Nx



Absolute discretion

Signia's tiny Pure® 10 sits nearly invisibly behind your ear. For a maximum of convenience and discretion, it can be remotely controlled via your smartphone.



Pure simplicity for hassle-free hearing

With no buttons required for manual adjustments, the fully automatic Pure 10 gives you the ultimate in simplicity and superior sound quality.

IP67 rated† – resistant to moisture, sweat, and dirt.

Hearing aids shown in actual size and may feature Color Exchange Kits.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The *Bluetooth®* word mark and logos are owned by the Bluetooth SIG, Inc., and any use of such marks by Sivantos Inc. is under license. Other trademarks and trade names are those of their respective owners.

Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. All other trademarks remain the property of their respective owners.

Apple App Store is a trademark of Apple Inc.

iPhone is a trademark of Apple Inc., registered in the U.S. and other countries.

iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

Android and Google Play are trademarks of Google Inc.

† Achieved IP67 rating per IEC 60529 standard. Water-resistant device can be completely submerged in water up to 3 feet for 30 minutes and dust will not interfere with the satisfactory operation of the device.

†† Achieved IP68 rating per IEC 60529 standard.

Legal Manufacturer

Sivantos, Inc.
10 Constitution Avenue
Piscataway, NJ 08854
Phone: (800) 766-4500
Fax: (732) 562-6696
www.sivantos.com

Copyright © 2018 Signia GmbH. All rights reserved.

signiausa.com SH201804 1.0 2/18 22B771 2K SI/18365-18



signia
Life sounds brilliant.

Replicating nature.

Signia Pure Nx hearing aids deliver the sound that nature intended.



signiausa.com/SigniaNx



myHearing App.

Your myHearing App

Your hearing care professional

Signia Pure Nx hearing aids.

Your natural choice for performance and hearing.

Whether you're in a meeting, taking an important business call, preparing a presentation, or meeting a friend for lunch, with Signia's new hearing aids you can enjoy our most natural sound experience along with superior connectivity and audio streaming from your smartphone and TV.

Hear the sound of your own voice as nature intended.

Reading a book to a child, having a little chat with your neighbors, looking at holiday photos with friends - Signia's hearing aids enable you to fully enjoy all of life's wonderful moments.

For the first time you can enjoy our most natural hearing experience and our highest hearing performance to perceive your own voice in perfect balance with your environment.

Please note: Features may vary depending on the configuration of the individual hearing aid. Your hearing care professional will be happy to assist you.

Enjoy direct streaming of phone calls, music, and TV.

Stream audio from a telephone or video conference directly into your ears, discreetly adjust the volume in a meeting, or listen to your navigation app for the fastest route to your destination - all you need is your smartphone.

Using the myControl™ App*, you can stream audio from various sources directly into your Pure 312 and Pure 13 hearing aids and control them remotely.

The touchControl™ App lets you change Pure 10's settings remotely and personalize your hearing aids to suit your individual preferences.

Further facts and figures.



Signia tinnitus solutions
Tune out tinnitus with the unique built-in tinnitus solutions from Signia. Ask your hearing care professional which solution is best for you.



Optional Telecoil
The battery door of Signia's Pure 13 Nx can be equipped with a Telecoil. In places with an induction loop, this helps you understand speakers more clearly. Simply select the T-program wherever you see the sign on the left.



Colors
Signia Pure Nx hearing aids come in a range of natural colors that are designed to match a variety of hair and skin tones. Individualize your hearing aids by choosing from a range of elegant options. Available colors vary depending on the model.

Personal care, anywhere.

Your journey back to effortless hearing is even easier with the myHearing™ App* for your smartphone. It lets you familiarize yourself with your new hearing aids in as little time as possible. You can follow special listening exercises, consult the user guide, and connect to your hearing care professional on the go.



* Compatible with Android and iOS. The apps are free of charge.

Find out more on myhearing.com or ask your hearing care professional for more information on Signia Nx hearing aids, useful accessories, and the myHearing App.