



ANIMAL HEALTH TEST OPTIONS

<p>Bartonella ePCR™ Single Blood or Single Fluid Draw</p> <p>Enrichment culture in patented BAPGM media provides enhanced detection of Bartonella DNA missed by standard PCR detection methods. Genus-level primers detect the presence of Bartonella DNA occurring at a detectable level. All positive results are verified by sequencing.</p>	<p>Specimen Requirements:</p> <p>Storage Instructions:</p> <p>Results Provided:</p> <p>Cost:</p>	<p>2-4 ml blood in purple top tube (EDTA)</p> <p>Refrigerated at 2-8C or frozen at minus 20C</p> <p>2-3 weeks</p> <p>\$125</p>
<p>Bartonella ePCR™ Triple Blood Draw</p> <p>Enrichment culture in patented BAPGM media provides enhanced detection of Bartonella DNA missed by standard PCR detection methods. Genus-level primers detect the presence of Bartonella DNA occurring at a detectable level. All positive results are verified by sequencing.</p>	<p>Specimen Requirements:</p> <p>Storage Instructions:</p> <p>Results Provided:</p> <p>Cost:</p>	<p>3 draws over 1 week (M-W-F) of: 4 ml blood in purple top tube (EDTA)</p> <p>Refrigerated at 2-8C or frozen at minus 20C</p> <p>2-3 weeks</p> <p>\$250</p>
<p>Bartonella ePCR™ Fresh/Frozen Tissue</p> <p>Pre-culture PCR of specimen followed by 10-day BAPGM enrichment culture and post-enrichment PCR. Sequence verification of any positive PCR results.</p>	<p>Specimen Requirements:</p> <p>Storage Instructions:</p> <p>Results Provided:</p> <p>Cost:</p>	<p>0.5 to 1 gram in red top or cryotube hydrated with sterile saline</p> <p>Refrigerated at 2-8C or frozen at minus 20C</p> <p>2-3 weeks</p> <p>\$200</p>
<p>Bartonella PCR Paraffin-Embedded Tissue</p> <p>Standard PCR is available for testing samples that cannot be cultured. Genus-level primers detect the presence of Bartonella DNA occurring at a detectable level. DNA sequencing and identification are included for positive results.</p>	<p>Specimen Requirements:</p> <p>Storage Instructions:</p> <p>Results Provided:</p> <p>CPT Coding:</p> <p>Cost:</p>	<p>2 cuts of minimum 50 microns thickness each, in Eppendorf tube</p> <p>Ambient</p> <p>1 week</p> <p>87471x1</p> <p>\$125</p>
<p>Bartonella Serology Panel</p> <p>Tests for antibodies against <i>Bartonella henselae</i>, <i>Bartonella vinsonii subsp. berkhoffii</i> and <i>Bartonella koehlerae</i>. Requires 3ml serum. Results reported within 1 week.</p>	<p>Specimen Requirements:</p> <p>Storage Instructions:</p> <p>Results Provided:</p> <p>Cost:</p>	<p>2ml-3ml of serum in a red top tube or in a tiger top tube</p> <p>Refrigerated at 2-8C or frozen at minus 20C</p> <p>1 week</p> <p>\$60</p>