

Which Stethoscope Software for Your Telemedicine System?

For RNK Products PCP Line of Stethoscopes	Auscultation Anywhere w/sSOIP Software	PCP-SSP Software
Stethoscope Compatibility	PCP-1, PCP-USB	PCP-USB
FDA Cleared for Telemedicine – Class II Medical Device	Yes	Yes
Separate Communications Channel for Stethoscope Signal	Yes	No
Uses Audio Channel of Video Conferencing System	No	Yes
Stethoscope Signal Can Remain Separate from Local Room Audio	Yes – By virtue of separate audio channel	Yes – By auto mute feature for local room audio
Stethoscope Signal Can Remain Separate from Remote Room Audio	Yes	No
Software Required at Local Station	Yes	Yes
Software Required at Remote (Physician) Station	Yes	No
PC Platform for Software	Windows OS	Windows OS
iPad, iPhone Usable by Physician	No	Possibly – If VC Software Supports IOS App
Local Monitoring of Stethoscope	Yes	Yes
Filters to Distinguish Heart and Lung Sounds	No Filter (20-2000 hz) 250 hz 500 hz 800 hz 1000 hz	No Filter (20-2000 hz) 250 hz 500 hz 800 hz 1000 hz Lungs (250 hz high pass)
Stethoscope Call Initiation	Remote (Physician) Site	Remote or Local Site Can Initiate VC Call
Recording Capability	Yes	No – Under Development
Physician Site License Fee	Yes	No
Patient Site License Fee	Yes	Yes
Annual Network Service Fee	Yes	No